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Ref. 1 ISO-20776. Clinical laboratory testing and in vitro diagnostic test systems - Susceptibility testing of infectious agents and evaluation of performance of antimicrobial susceptibility test devices - Part 2: Evaluation of performance of antimicrobial susceptibility test devices. International Organization for Standardization, 2007.

Ref. 2. U.S. Department of Health and Human Services Food and Drug Administration, Center for Devices and Radiological Health. Guidance for Industry and FDA. Class II Special Controls Guidance Document: Antimicrobial Susceptibility Test (AST) Systems, August 28, 2009.

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Clinical and Laboratory Standards Institute (CLSI). **2015** Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard - Tenth Edition. CLSI Document M07-A010, Wayne, PA.

CLSI-EUCAST Polymyxin Breakpoints Working Group. **2016**. Recommendations for MIC determination of colistin (polymyxin E) as recommended by the joint CLSI271 EUCAST Polymyxin Breakpoints Working Group. Available at:
http://www.eucast.org/fileadmin/src/media/PDFs/EUCAST_files/General_documents/Recommendations_for_MIC_determination_of_colistin_March_2016.pdf. Accessed 4 August **2017**

Clinical and Laboratory Standards Institute (CLSI). **2017**. Performance standards for antimicrobial testing. CLSI document M100-S27, Wayne, PA.

Humphries RM. Susceptibility Testing of the Polymyxins:Where Are We Now? **2015**. Pharmacotherapy. 35:22-27.

EUCAST. **2017**. Clinical Breakpoints version 7.1 (Update 2017-03-13). Available at:
http://www.eucast.org/clinical_breakpoints/. Accessed 4 August 2017.

Chew KL, La MV, Lin RTP, Teo JWP. 2017. Colistin and polymyxin B susceptibility testing for carbapenem-resistant and *mcr* -positive *Enterobacteriaceae*: Comparison of Sensititre, Microscan, Vitek 2, and Etest with broth microdilution. *J Clin Microbiol*.doi:10.1128/JCM.00268-17.

EUCAST. **2017**. Warnings concerning antimicrobial susceptibility testing products or procedures. Antimicrobial susceptibility testing of colistin - problems detected with several commercially available products. Available at:http://www.eucast.org/ast_of_bacteria/warnings/. Accessed 4 August 2017.

Rojas LJ, Salim M, Cober E, Richter SS, Perez F, Salata RA, Kalayjian RC, Watkins RR, Marshall S, Rudin SD, Domitrovic TN, Hujer AM, Hujer KM, Doi Y, Kaye KS, Evans S, Fowler VG Jr, Bonomo RA, van Duin D; Antibacterial Resistance Leadership Group. 2017. Colistin Resistance in Carbapenem-Resistant *Klebsiella pneumoniae*: Laboratory Detection and Impact on Mortality. *Clin Infect Dis*. **64:711-718**

Dafopoulou K, Zarkotou O, Dimitroulia E, Hadjichristodoulou C, Gennimata V, Pournaras S, Tsakris A. **2015**. Comparative Evaluation of Colistin Susceptibility Testing Methods among Carbapenem-Nonsusceptible *Klebsiella pneumoniae* and *Acinetobacter baumannii* Clinical Isolates. *Antimicrob Agents Chemother*. Aug;59(8):4625-30.

Vasoo S. 2017. Susceptibility Testing for the Polymyxins: Two Steps Back, Three Steps Forward?. *J Clin Microbiol*. Sep;55(9):2573-2582. doi: 10.1128/JCM.00888-17.

Jerke K., Lee M., Romney M. Humphries R. **2016**. Polymyxin Susceptibility Testing: a Cold Case Reopened. *Clinical Microbiology Newsletter* 38:69-77.

Vourli S, Dafopoulou K, Vrioni G, Tsakris A, Pournaras S. **2017**. Evaluation of two automated systems for colistin susceptibility testing of carbapenem-resistant *Acinetobacter baumannii* clinical isolates. *J Antimicrob Chemother.* Sep 1;72(9):2528-2530

Poirel L, Jayol A, Nordmann P. 2017. Polymyxins: Antibacterial Activity, Susceptibility Testing, and Resistance Mechanisms Encoded by Plasmids or Chromosomes. Clin Microbiol Rev. 2017;30:557-596.

Liu Y, Wang Y, Walsh TR, Yi LX, Zhang R, Spencer J, Doi Y, Tian G, Dong B, Huang X, Yu LF, Gu D, Ren H, Chen X, Lv L, He D, Zhou H, Liang Z, Liu JH, Shen J. **2016**. Emergence of plasmid-mediated colistin resistance mechanism MCR-1 in animals and human beings in China : a microbiological and molecular biological study. *Lancet Infect Dis.* 16:161-168.

Tijet N, Faccone D, Rapoport M, Seah C, Pasterán F, Ceriana P, Albornoz E, Corso A, Petroni A, Melano RG. 2017. Molecular characteristics of mcr-1-carrying plasmids and new mcr-1 variant recovered from polyclonal clinical Escherichia coli from Argentina and Canada. PLoS One. Jul 5;12(7):e0180347.

Insert for Kit for Detection of COLISTIN resistance in Enterobacteriaceae, *P.aeruginosa* and *Acinetobacter* (98018). <http://www.rosco.dk/> Accessed 12 June **2017**

Dominguez JE, Figueroa Espinosa RA, Redondo LM, Cejas D, Gutkind GO, Chacana PA, Di Conza JA, Fernández-Miyakawa ME. **2017**. Plasmid-mediated colistin resistance in *Escherichia coli* recovered from healthy poultry. *Rev Argent Microbiol.* 49(3):297-298.

Liakopoulos A, Mevius DJ, Olsen B, Bonnedahl J. **2016**. The colistin resistance mcr-1 gene is going wild. *J Antimicrob Chemother.* 1(8):2335-6. doi: 10.1093/jac/dkw262.

Rapoport M, Faccone D, Pasteran F, Ceriana P, Albornoz E, Petroni A; MCR Group, Corso A. **2016**. First Description of mcr-1-Mediated Colistin Resistance in Human Infections Caused by *Escherichia coli* in Latin America. *Antimicrob Agents Chemother* 60(7):4412-3.

Chew KL, La MV, Lin RTP, Teo JWP. 2017. Colistin and Polymyxin B Susceptibility Testing for Carbapenem-Resistant and mcr-Positive Enterobacteriaceae: Comparison of Sensititre, MicroScan, Vitek 2, and Etest with Broth Microdilution. J Clin Microbiol. Sep;55(9):2609-2616.

Hindler J, Humphries RM. **2013**. Colistin MIC variability by method for contemporary clinical isolates of multidrug-resistant Gram-negative bacilli. *Clin Microbiol.* Jun;51(6):1678-84

Giamarellou H. **2016**. Epidemiology of infections caused by polymyxin-resistant pathogens. *Int J Antimicrob Agents.* 48:614-621

Nordmann P, Jayol A., Poirel L. Rapid Detection of Polymyxin Resistance in Enterobacteriaceae **2016**. *Emerging Infectious Diseases.* 22: 1038-43.

Jayol A, Nordmann P, Lehours P, Poirel L, Dubois V. **2017**. Comparison of methods for detection of plasmid-mediated and chromosomally encoded colistin resistance in Enterobacteriaceae. Clin Microbiol Infect. 17:30291-4. doi: 10.1016.

Yin W, Li H, Shen Y, Liu Z, Wang S, Shen Z, Zhang R, Walsh TR, Shen J1, Wang Y. **2017**. Novel Plasmid-Mediated Colistin Resistance Gene *mcr-3* in *Escherichia coli*. MBio. Jun 27;8(3). pii: e00543-17. doi: 10.1128/mBio.00543-17.

Al-Tawfiq J., Laxminarayan R, Mendelson M. **2016**. How should we respond to the emergence of plasmid-mediated colistin-resistance in humans and animals? IJID. doi:10.1016

Jeannot K, Bolard A., Plesiat P. Resistance to polymyxins in Gram-negative organisms. **2017**. Int J Antimicrob Agents.;49:526-535.

ASM Microbe 2017, 01 al 05 de Junio, New Orleans, EEUU

- **Epidemiology of *mcr-1*-producing Enterobacteriaceae clinical isolates from Argentina. Faccone D., M. Rapoport, E. Albornoz, P. Ceriana, C. Lucero, P. Gagetti, A. Petroni, MCR-Group, F. Pasteran, A. Corso. S-135/2799**

- **Screening of Emerging Resistant Mechanisms (KPC, NDM, OXA-48-Like and *mcr*) in Enterobacteriaceae with MicroScan Automated System. F. Pasteran, P. Ceriana, M. Rapoport, E. Albornoz, C. Lucero, P. Gagetti, G. Baldoni, A. Corso. F- 421/3734.**

- **Rapid Polymyxin NP Test for Rapid Detection of Polymyxin-Resistant *Enterobacteriaceae*. N. Kieffer, L. Poirel, P. Nordmann. S- 68/2117**

- **Underestimation of Colistin Resistance among *Acinetobacter baumannii* Clinical Isolates by Two Automated Systems Compared with Agar and Microbroth Dilution Methods .S. Vourli, M. Beredaki, K. Dafopoulou, J. Meletiadiis, A. Vasilakopoulou, A. Tsakris, S. Pournaras. S- 144/4500**

- **Comparison of Agar Dilution and Broth Microdilution for Colistin with a Diverse Set of Gram-Negative Isolates. F. Basit, C. Shubert, S. Messina-Powell, M. F- CPHM LB1/6602**

- **Comparison of Etest and Broth Microdilution Methods for Detection of Colistin Resistance in Clinical Isolates. S. Simar, D. Sibley, D. Ashcraft, G. Pankey; F- 398/5341**

- **Testing of USA *mcr-1* Colistin-Resistant *E. coli* on MicroScan Panels. B. Zimmer, E. Fernandez, F. Popoola, R. Kwong; FRIDAY - 406/ 3478**

- **Evaluation of a Broth Microdilution Plate for Determining the Susceptibility of Carbapenem-Resistant *Klebsiella pneumoniae* to Polymyxins . J. Karichu¹, J. Otiso¹, H. Van Heule¹, G. Keller¹, E. Cober¹, L. J. Rojas², A. M. Hujer², K. M. Hujer², S. Marshall², F. Perez², S. D. Rudin², T. N. Domitrovic², K. S. Kaye³, R. S. Salata⁴, D. van Duin⁵, R. A. Bonomo², S. S. Richter. F- 432/4643**

- **Comparison of Different Commercial Methods for Polymyxin Susceptibility Testing. P. C. M. Koga, A. M. Doi, K. Santiago, T. B. Melo, J. Pasternak, M. D. V. Martino. F - 430/ 4737**