

## Colistin Drop Test

Version

EN 1.0

Last update: May 2020

Pages: 3



Strengths	<ul style="list-style-type: none"> <li>Colistin solution is stable for more than 1 year in the refrigerator</li> <li>Test can be carried out on the same plate on which disks have been placed to carry-out an antibiogram</li> <li>Is capable to detect colistin resistance mediated by traditional (chromosomic) as well as transferable plasmidic (<i>mcr</i>) mechanisms.</li> </ul>
Limitations	<ul style="list-style-type: none"> <li>Cation-adjusted Mueller-Hinton Broth (CA-MHB) required</li> <li>Validated for Difco®, Oxoid® and BD Phoenix™ AST broth (Cat. 246003) Mueller Hinton Broth</li> <li>Validated for Oxoid® and BBL® paper disks</li> </ul>
Organism group	<i>Enterobacterales</i> , <i>Pseudomonas aeruginosa</i> and <i>Acinetobacter</i> spp.
Medium	<ul style="list-style-type: none"> <li>Agar Mueller-Hinton plates (as used for the disk diffusion method)</li> <li>CA-MHB</li> </ul>
Antimicrobial concentration	Colistin sulfate 16 µg/ml
Source of antimicrobial	10-µg colistin paper disks <i>Alternatively, colistin sulphate can be used</i>
Inoculum	<ul style="list-style-type: none"> <li>Using a loop or a swab, pick 3-5 colonies from a fresh (&lt;24 h) culture (selective or non-selective media can be used) and transfer to sterile saline or broth.</li> <li>Adjust turbidity to equivalent of a 0.5 Mc Farland turbidity standard</li> </ul>
Test procedure	<p>A) <i>Preparation of the colistin sulfate 16 µg/ml solution:</i></p> <ul style="list-style-type: none"> <li>Add 8 colistin disks (10-µg) in a glass tube containing 5 ml of CA-MHB.</li> <li>Incubated at room temperature for at least 30 minutes but not longer than 60 minutes to allow colistin to elute from the disks.</li> <li>Remove aseptically the disks and discard.</li> <li>The remaining broth is the working solution (colistin final concentration of 16 µg/ml).</li> <li>Colistin solution can be stored for up to 1 year at 4 °C.</li> <li>It is suggested to monthly control the sterility of the solution by overnight incubation at 35 ° C of 100 µl</li> </ul>

	<p>and observation of turbidity. If this happens, discard the solution. Alternatively, it can be re-sterilized by filtration (0.2µm filter pore size required).</p> <p><b>B) Drop Test procedure:</b></p> <ul style="list-style-type: none"> <li>• Swab the surface plate with a 0.5 Mc Farland suspension of the bacterial isolate.</li> <li>• Preferably, leave the plate to rest with its lid on for approx. 15 min to dry the inoculum.</li> <li>• Using a pipette or a calibrated loop, add a 10 µl drop of colistin sulfate solution (16 µg/ml) over the plate surface already inoculated.</li> <li>• If necessary, reserve the plates for 15 minutes at room temperature to allow complete absorption of the drop, then invert the plate.</li> <li>• Incubate for 16 to 18 hours at 35 °C.</li> <li>• After incubation, examine the purity plate to ensure inoculum was pure.</li> <li>• Observe the presence or absence of an inhibition zone.</li> </ul>
Results	<ul style="list-style-type: none"> <li>• COLISTIN-SUSCEPTIBLE: Clear inhibition zone, of any diameter.</li> <li>• COLISTIN-RESISTANT: <ul style="list-style-type: none"> <li>○ Absence of inhibition zone or</li> <li>○ Defined colonies or haze growth within the inhibition zone (indicative of hetero-resistant subpopulations)</li> </ul> </li> </ul> <p><i>See figure in Annex for clarification</i></p>
Additional testing and reporting	
QC recommendations	<ul style="list-style-type: none"> <li>• Colistin susceptible: <i>Escherichia coli</i> ATCC 25922</li> <li>• Colistin resistant: <i>E. coli mcr</i> producer</li> <li>• <i>The quality of each batch of Agar Mueller-Hinton plates and CA-MHB should be previously controlled according to the latest edition of CLSI<sup>1,2</sup> or EUCAST<sup>3</sup></i></li> </ul>
Reference	<p><sup>1</sup> Clinical and Laboratory Standards Institute. M02: Performance Standards for Antimicrobial Disk Susceptibility Tests.</p> <p><sup>2</sup> Clinical and Laboratory Standards Institute M07: Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically</p> <p><sup>3</sup> European Committee on Antimicrobial Susceptibility Testing. Media preparation for EUCAST disk diffusion testing and for determination of MIC values by the broth microdilution method</p> <p><sup>4</sup> Development and validation of simple tests (agar spot, colistin drop, 1ml-broth disk elution MIC and tablet pre-diffusion) as an alternative to improve accuracy in screening chromosomal and plasmid-mediated colistin resistance in GNB. F. Pasteran, D. Danze, C. Cabrera, C. Lucero, A. Menocal, E. Albornoz, I. Castillo, M. Rapoport, P. Ceriana, P. Gagetti, A. Corso. O0952. 28<sup>o</sup> ECCMID, 2018</p>

Annex.

**Figure. Example of colistin susceptible and colistin resistant bacteria.**

